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The Effect on HPLC Fingerprints of Flos Genkwa Prepared by Macroporous Resin

LI Rao-rao^{1*}, WANG Cai-fang^{1,2}, HUANG Lan-lan^{1,3}, ZHONG Li-qun⁴, YUAN Si-tong¹

(1. Institute of Chinese Materia Medica, China Academy of Chinese Medical Sciences, Beijing 100700, China;
2. College of Pharmacy, Zhengzhou University, Zhengzhou 450052, China;
3. Jiangxi College of Traditional Chinese Medicine, Nanchang 330006, China;
4. Dongzhimen Hospital, Beijing University of Traditional Chinese Medicine, Beijing 100700, China)

[Abstract] **Objective:** To analyze the effect on HPLC fingerprint of Flos Genkwa which samples prepared by macroporous resin. **Method:** The fingerprints of 17 batches of Flos Genkwa were determined. The samples were extracted with methanol, then solvent were distilled off under reduced pressure. The residue were solved by water. Macroporous resin column chromatography was used for separating the water solution. Three parts of eluates, 20%, 50% and 95% ethanol eluate were obtained, separately. Then they were analyzed by high performance liquid chromatography. Chromatographic conditions: chromatographic column was Diamonsil[®] C₁₈ (4.6 mm × 250 mm, 5 μm), mobile phase was methanol-water in gradient elute, flow rate was 1 mL·min⁻¹, column temperature was 35 ℃, detective wavelength was set at 332 nm, record time was 40 min. **Results:** 7 common peaks were obtained in methanol extract, while 11 common peaks were displayed after prepared by macroporous resin. **Conclusion:** It can show more information in HPLC fingerprint which samples prepared by macroporous resin. This is a useful method in sample preparation.

[Key words] Flos Genkwa; macroporous resin preparation method; HPLC; fingerprint

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