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Determination of Isoniazid and Pyrazinamide in Human Plasma by HPLC

Zhao Guanren, Ma Jun, Feng Duanhao (Department of Pharmacy, the 309th Hospital of PLA, Beijing 100091, China)

ABSTRACT **Objective:** To determine the concentration of isoniazid and pyrazinamide in human plasma by RP-HPLC. **Method:** The HPLC method was carried out on a Diamonsil C₁₈ (250 mm × 4.6 mm, 5 μm) column with 0.05 mol · L⁻¹ KH₂PO₄-acetonitrile (89: 11, adjusting pH to 6.0 with triethylamine) as the mobile phase. The flow rate was 0.9 ml · min⁻¹, and the column temperature was ambient with the detection wavelength of 254 nm. **Result:** The linear range of isoniazid was 0.5-30 μg · ml⁻¹ ($r=0.9984$), and that of pyrazinamide was 1-50 μg · ml⁻¹ ($r=0.9981$). The extraction recovery was 90.33% and 103.01%, respectively. **Conclusion:** The method is sensitive, simple and rapid. It can be used in the determination of isoniazid and pyrazinamide in human plasma. **KEY WORDS** HPLC; Isoniazid; Pyrazinamide; Plasma concentration

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Determination of Fluorouracil and Related Substances in Compound Fluorouracil Injections

Zhang Di¹, Guo Hongwei¹, Sun Kuan¹, Sun Lingling¹, Yang Cuo², Zhao Chunjie² (1. Liaoning Institute for Food and Drug Control, Shenyang 110023, China; 2. Shenyang Pharmaceutical University)

ABSTRACT Objective: To establish a RP-HPLC method for the determination of fluorouracil and its related substances in compound fluorouracil injections. **Method:** The determination was performed by HPLC with C₁₈ column. The mobile phase was water (adjusting pH to 3.5 with 0.05 mol·L⁻¹ phosphoric acid) - methanol (95:5), and the UV detection wavelength was 265 nm. **Result:** Fluorouracil could be completely separated with the other impurities. The linear rang was 0.038 92 ~ 1.946 μg ($r = 0.999\ 9$, $n = 6$), the average recovery was 99.2% (RSD = 0.5%) with the detection limit of 0.8 ng. **Conclusion:** This method is selective and accurate, and suitable for the quality control of compound fluorouracil injections.

KEY WORDS Compound fluorouracil injections; Content determination; Related substances

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